

# Cable Management: Mainstream CAPNOSTAT 5 Sensor


*The following instructions demonstrate proper cable management for the CAPNOSTAT 5 mainstream sensor.*



Always coil the cable in the direction of the natural cable coil, greater than 3 inches [7.6 cm] in diameter.



DO NOT kink the sensor cable, or bend the sensor cable less than 3 inches [7.6 cm] in diameter.



The use of cable ties is acceptable, and will aid in maintaining the diameter of the cable during storage. This illustration shows the positioning of two cable ties to avoid kinking the cable at the connector and sensor strain reliefs.



DO NOT coil the sensor cable less than 3 inches [7.6 cm] in diameter.




This illustration shows positioning of one cable tie. Position the cable tie to avoid sharp kinks at the connector and sensor strain reliefs.



DO NOT coil and band the cable, as this will kink the cable and bend the cable smaller than 3 inches [7.6 cm] in diameter.




DO NOT wrap the sensor cable around the airway adapter channel on the sensor.



DO NOT coil the sensor cable around a ventilator circuit support arm, monitor stand, or any object less than 3 inches [7.6 cm] in diameter.




DO NOT wrap the sensor cable around the sensor head.



DO NOT pinch and straighten the cable during cleaning, handling, or use.

# Cable Management: Sidestream LoFlo Sensor

*The following instructions demonstrate proper cable management for the sidestream LoFlo sensor.*



Always coil the cable in the direction of the natural cable coil, greater than 3 inches [7.6 cm] in diameter. The use of cable ties is acceptable, and will aid in maintaining the diameter of the cable during storage. Position the cable tie to avoid sharp kinks at the connector and sensor strain reliefs.



DO NOT wrap the cable around the sensor in any manner.



DO NOT kink the sensor cable, or bend the sensor cable less than 3 inches [7.6 cm] in diameter.



DO NOT wrap the cable around the optional pole mount in any manner.



DO NOT coil and band the cable, as this will kink the cable and bend the cable smaller than 3 inches [7.6 cm] in diameter.



DO NOT wrap the cable around the optional pole mount in any manner.




DO NOT fold the cable over on itself.



DO NOT coil the sensor cable around a ventilator circuit support arm, monitor stand, or any object less than 3 inches [7.6 cm] in diameter.



DO NOT wrap the cable around the sensor in any manner.



DO NOT pinch and straighten the cable during cleaning, handling, or use.