

veterinary capnography accessories




Mainstream and Multigas



IRMA CO₂ analyzer 4420505
IRMA AX+ multigas analyzer 4420506

Cable Length: 2.5m



IRMA Airway Adapter Single Patient Use

Large	4420511	Pack of 25
Small	4420521	Pack of 10

Sidestream

Short - Long Term (Single Patient Use)



Sampling Line
4420550 Box of 40
Length: 2m
Connections: Male Luer
Lock on both ends



Airway Adapter (Large/Small) 4420531 Box of 25

Alternative 15mm Airway Adapters	Boxes of 10	
I.D. 2.0mm 4420570	I.D. 4.0mm 4420574	I.D. 8.0mm 4420578
I.D. 2.5mm 4420571	I.D. 5.0mm 4420575	I.D. 9.0mm 4420579
I.D. 3.0mm 4420572	I.D. 6.0mm 4420576	I.D. 10.0mm 4420580
I.D. 3.5mm 4420573	I.D. 7.0mm 4420577	



**Nomo Adapter
(Long Term)**
4420525
Length: 22cm

with female
Luer lock
connector



Length: 2m

Nasal Sampling Cannula

Large	4420710	Box of 25
Small	4420711	Box of 25

Short Term/Emergency (Single Patient Use)



Length: 2m

Nasal Sampling Line

For triage and emergency applications

4420560 Box of 40



Thames Medical Ltd.
Tel: +44 (0)1903 522911
www.thamesmedical.com