

## Your Lutech Datalys Multiparameter Monitor Disposable item Replacement Schedule

- Cleaning and Disinfecting:** After every use:  
Use only **non-chlorine** based cleaner (eg: TECCare Control)  
**Never** immerse items in liquid
- Airway Adaptors Large:** Replace at least every 1 Month
- Airway Adaptors Small:** Replace at least every 1 Month
- Monitor:** Service & Calibration every 12 Months

### Replacement Schedule

Month	Airway Adaptor - Large	Airway Adaptor - Small	Service & Calibration
January			
February			
March			
April			
May			
June			
July			
August			
September			
October			
November			
December			

**Contact Thames Medical to arrange your Annual Service & Calibration  
Annual Disposable and Cleaning Packs available**

Thames Medical recommends:

