

Your VM-2500 Mainstream Capnograph

Disposable item Replacement Schedule

Cleaning and Disinfecting:	After every use: Use only non-chlorine based cleaner (eg: TECCare Control) Never immerse items in liquid
Airway Adaptors:	Replace at least every 1 Month
Monitor:	Monitor Zeroed & Calibration every 12 Months

Replacement Schedule

Month	Irma Airway Adaptor Adult/Paediatric (large)	Irma Airway Adaptor Infant (small)	Monitor Zeroed/ Calibration
January			
February			
March			
April			
May			
June			
July			
August			
September			
October			
November			
December			

Contact Thames Medical to arrange your Service & Calibration
Disposable and Cleaning Packs available

Doc Ref: TML-QC-P-CAPS-M-V2.0